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# FreezeBox Instructions for use

## Overview

Ice-free Freeze box is a new type of energy-saving, environmentally friendly, safer, and standardized with no additives, cooling equipment. It uses no electricity, no ice and brings low temperature protection for biological samples.

Ice core aluminum alloy shell with high thermal conductivity of the alloy module to ensure that the sample can be quickly cooled, the temperature difference between the sample holes less than 0.1°C, it can maintain a uniform temperature between the samples.

The Freeze box is Ice-free, lightweight, easy to carry, ideal for collection of biological samples, it also can be used in helping to preserve the unstable nucleic acid samples. Its compact design and ice-free technology make it ideal for tissue culture incubator, bio-safety cabinet and any other application with limited space.

The freeze box can choose with different ice core, and a variety of thermal test tube modules for usage in different samples and temperatures.

## Product Features

1. Can bring low temperature protection for biological samples all mes;
2. The uniformity of the temperature between the holes can ensure consistency of the sample test results;
3. Optional different ice cores selection for different storage temperatures;
4. Applicable to a variety of thermal pipes and modules;
5. Easy to clean and disinfect;
6. Lightweight, easy to carry, suitable for the collection of biological samples;
7. Highly durable, suitable for frequent and repeated use.

## Use method

By selecting from a number of ice core and cooling options, samples can be maintained within a variety of temperature bands that range from -78°C to 4°C.

Temperature Range	Cooling Sources	Insulation Time (Open cover)	Insulation Time (Close cover)
0.5°C~4°C	Ice core CL	Over 4 hours	Over 10 hours
-18°C~-4°C	Ice core CF	Over 3 hours	Over 6 hours
-30°C~-20°C	Dry ice with heat shield	Over 5 hours	Over 6 hours
-78°C	Dry ice	Over 5 hours	Over 6 hours


NOTE: Use of the heat shield (open cover) to protect exposed surface of Tube module from environmental heat influx will provide optimal cooling duration. The heat shield is the choice and buy.

Actual performance may vary depending upon ambient temperature, start-up parameters, sample load, initial sample temperature, air currents, radiant energy sources and other conditions.

### 1. 0.5°C~4°C (Ice core CL)

Freeze ice core CL in a -20°C freezer for a minimum of four hours. Pre-chill the tube module to 0°C to 4°C in a refrigerator for approximately 30 minutes. Remove ice core from freezer and allow ice core to sit at room temperature approximately 10 minutes. Place a frozen ice core in the freezebox base followed by a tube module. The module temperature should equilibrate between 0.5°C to 4°C. Keep the lid on the freezebox to prolong the cooling duration; use heat shield to extend ice core service interval or conditions of higher environmental temperatures.


Note: Ice core and tube module of the processing can also didn't go to other way, as shown in the figure below

### 2. -18°C~-4°C (Ice core CL)

Freeze ice core CF in a -20°C freezer for a minimum of 6 hours. Pre-chill the tube module in the -20 freezer for approximately one-half hour. Place a frozen ice core in the freezebox followed by the pre-chilled tube module. Keep the lid on the freezebox to prolong the cooling duration; use heat shield to extend ice core service interval or conditions of higher environmental temperatures.


Note: Ice core and tube module of the processing can also didn't go to other way, as shown in the figure below

### 3. -78°C~-20°C (Dry ice)

Fill the freezebox cavity with pulverized dry ice level with the bottom of the finger grip recesses as shown. Place the heat shield onto the dry ice. Place the tube module on top of the heat shield. Tube module temperatures will quickly equilibrate to approximately -20°C and solely drop to -30°C as the dry ice sublimates. Place the tube module directly onto dry ice. Tube module temperatures will quickly reach dry ice temperatures in approximately 8 minutes. Closing the lid will not further decrease the tube module temperatures, however it will extend the cooling duration.

Note: Specific data, as shown in the figure below

## 1. Care and cleaning

Freezebox is constructed from high density closed cell polyethylene foam. Freezebox is compatible with prolonged cryogenic temperature exposure. The foam may be cleaned by water and mild soap. Rinse thoroughly. The freezebox is resistant to alcohols and 10% bleach solutions. Avoid abrasive or sharp objects. Do not use the freezebox for pulverizing dry ice. Do not autoclave. Maximum temperature exposure: 60. Avoid exposure to UV light sources. The freezebox is not recommended for shipping of temperature-sensitive samples.

## 2. Test Equipment

Temperature probe with small thermal mass such as a thermocouple probe  
Digital thermometer to match probe.  
Representative closed sample tube with cap.

## 3. Test procedure

Drill a small hole in the center of the sample tube cap with a diameter that provides a snug fit for the temperature probe wire so there will be minimal ambient air influx into the tube

Fill the sample tube with 1 ml of water

Place the cap on sample tube and insert temperature probe

Place the instrumented sample tube in any well of the tube module

Rack temperatures can be measured by direct contact of the probe in any of the rack wells

Caution: The product described here are intended for the exclusive use of the methods described. Require the use of dry ice. Direct skin contact with dry ice or metal components that have been touching dry ice can cause freezing injury. Always use appropriate protective equipment for eyes and skin when handling dry ice and cold metal components.