
NEST®

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
NEST®


NEST Multi-Layer Cell Culture Flask Culturing Technique

NEST multi-layer flasks provide 174 cm² cell growth surface area per layer. The tissue culture treated NEST multi-layer flasks enable you to grow more cells faster and easier, thereby making your cell culture workflow more productive, easier, and in a more convenient manner.


ADDING MEDIA AND PREPARING CELL SUSPENSION

1. Mix cell suspension with medium: Prepare cell suspension of required concentration in a container and mix. Recommended volume is about 30-50 mL per layer.
2. Add the mixed liquid into the multi-layer flask slowly with serological pipettes. To avoid foams and bubbles, allow liquid stream to flow along the slope of the multi-layer flask. (Save a little liquid in pipette for each dosing to prevent bubbles.)
Tips: A 10 mL pipette allows media to be dispensed at the bottom of the vessel. A 25 mL pipette allows media to be dispensed at the spot just past the NEST logo.


3. Hold the multi-layer flask upright with the logo facing you and tilt clockwise to a 45° angle on a flat work surface to partition the liquid into each layer.
4. While holding the multi-layer flask at a 45° angle, gently lay it flat onto the work surface with the logo facing up.


5. After placing the multi-layer flask flat on a work surface, gently rock back and forth and side-to-side to distribute cells evenly onto culture surfaces.
Tips: Take care to avoid foaming of medium, and not to spill liquid from each layer.
6. Gently horizontally transfer and place the flask to a CO₂ incubator.

MEDIA REMOVAL

You may choose to either aspirate or pour the media from the multi-layer flask.

7. Aspirating Method:
To aspirate or remove media, with the NEST logo facing you, tilt the multi-layer flask counter-clockwise to a 45° angle while inverting the multi-layer flask toward you. Then, tilt multi-layer flask to the right, continuing to aspirate all residual media.
Tips: Aspirating media completely requires the usage of a NEST 2mL or 10mL serological pipette.

8. Pouring Method:
With the logo facing you, pour the media from the multi-layer flask.

CELL HARVESTING

9. Wash with buffer once and add dissociating reagent (≥ 5 mL per layer). Then, follow Steps 3-4 to distribute the dissociating reagent to each layer.
10. Neutralize with inactivating solution and mix it following Steps 3-4. Gently swirl to dislodge cells completely.
11. a). Follow Step 7 "Aspirating Method" protocol and collect cell suspension using a NEST 10mL serological pipette.
b). Follow Step 8 "Pouring Method". Pour the detached cell suspension into a NEST conical tube.
12. Rinse with additional wash buffer as needed.

Cat. No.	Description	/pack	/case
731001	5-layer cell culture flask, straight neck, plug seal cap, TC treated	1	8
731002	5-layer cell culture flask, straight neck, vent cap, TC treated	1	8
731301	3-layer cell culture flask, straight neck, plug seal cap, TC treated	1	12
731302	3-layer cell culture flask, straight neck, vent cap, TC treated	1	12

Wuxi NEST Biotechnology Co., Ltd.

Add: No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China, 214112
Tel: (86)510-6800 6788
E-mail: info@nest-wuxi.com
http://www.cell-nest.com

Nest Scientific Inc.

Add: Nest Scientific Inc., 3 Convery Blvd, Unit 600, Woodbridge NJ 07059, USA
Tel: (1)732-381-0268
E-mail: info@nestscientificusa.com
http://www.nestscientificusa.com

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