

移液器品牌	移液器型号	通道	通道规格 (μL)	301006	305006/ 302106	304006/ 303206	301016/ 301081	305016/302116/ 302081/305081	304016/ 303216/30308 1	
Axygen	Axypet	单通道	0.1-2	•			•			
			0.5-10	•			•			
			2-20		•			•		
			5-50		•			•		
			10-100		•			•		
			20-200		•			•		
			100-1000			•			•	
		多通道	1-10	•				•		
			5-50		•			•		
			20-200		•			•		
50-300			•			•				
Biohit	Proline	单通道	0.1-2.5	•			•			
			0.2-10	•			•			
			0.5-10	•			•			
			5-50		•			•		
			5-100		•			•		
			10-250		•			•		
			10-500			•			•	
			20-200		•			•		
			50-1000			•			•	
			100-1000			•			•	
	多通道	0.2-10	•				•			
		5-100		•			•			
		25-250		•			•			
		50-250		•			•			
		5-100		•			•			
	e-Line	单通道	0.1-10	•				•		
			0.2-10	•				•		
			5-120		•				•	
			10-300		•				•	
			50-1000			•			•	
		多通道	0.2-10	•				•		
			5-120		•			•		
			10-300		•			•		
			5-100		•			•		
			30-300		•			•		
m-Line	单通道	0.1-3	•				•			
		0.5-10	•				•			
		10-100		•				•		
		20-100		•				•		
	多通道	0.5-10	•				•			
		5-100		•			•			
		30-300		•			•			
		100-1000			•			•		
BrandTech	Transferpet te	单通道	0.1-1	•			•			
			0.5-10	•			•			
			2-20		•			•		
			5-50		•			•		
			10-100		•			•		
			20-200		•			•		
			25-250		•			•		
			100-1000			•			•	
		多通道	0.5-10	•				•		
			2-20		•			•		
			2.5-25		•			•		
			5-50		•			•		
			10-100		•			•		
			20-200		•			•		
			30-300		•			•		

Capp Aero		单/多通道	0.2-2	•			•			
			0.5-10	•			•			
			5-50			•		•		
			10-100			•		•		
			25-200			•		•		
			50-300			•		•		
			100-1000				•		•	
Labnet	Labpette Discovery or Biopette™ Plus	单通道	0.1-2	•			•			
			0.5-10	•			•			
			2-20			•		•		
			10-100			•		•		
			20-200			•		•		
			100-1000				•		•	
	多通道	1-10	•				•			
		5-50			•		•			
		20-200			•		•			
		50-300			•		•			
	BioPette E	单通道	0.5-10	•			•			
			2-20	•			•			
			10-200			•		•		
			100-1000				•		•	
多通道		0.5-10	•				•			
		2-20	•				•			
		10-200			•		•			
		100-1200				•		•		
Matrix	Impact/Impact ²	单/多通道	2-125			•		•		
Nichryo/Oxford	Benchmate	单通道	0.5-10	•			•			
			2-20			•		•		
			10-50			•		•		
Rainin	EDP3	单通道	E-MIC-10,0.5-10	•			•			
			E2-25,2.5-25			•		•		
			E2-100,10-100			•		•		
			E2-250,25-250			•		•		
			E2-1000,100-1000				•		•	
	多通道	EP-M8-10,0.5-10	•				•			
		EP-M8-250,10-250			•		•			
	Latch Mode	单通道	R2,0.1-2	•			•			
			R10,0.2-10	•			•			
			R20,2-20			•		•		
			R100,10-100			•		•		
R200,20-200					•		•			
R1000,100-1000						•		•		
Socorex	Callbra/Accura	单通道	0.1-2	•			•			
			0.2-2	•			•			
			0.5-10	•			•			
			1-10	•			•			
			2-20			•		•		
			10-100			•		•		
			20-200			•		•		
			100-1000				•		•	
		多通道	1-10	•				•		
			5-50			•		•		
			10-100			•		•		
VWR	VWR and VWR Signature	单通道	0.1-2	•			•			
			0.5-10	•			•			
			2-20			•		•		
			5-50			•		•		
			10-100			•		•		
			20-200			•		•		
			50-250			•		•		

	Signature	多通道	100-1000			•			•	
			1-10	•			•			
			5-50		•			•		
			20-200		•			•		
			50-300		•			•		
Eppendorf	Electronic	单/多通道	0.5-10	•			•			
			5-100		•			•		
			2-200		•			•		
			50-1000			•			•	
	Reference	单/多通道	0.1-2.5	•			•			
			0.1-10	•			•			
			0.5-10	•			•			
			2-20		•			•		
			10-100		•			•		
			20-200		•			•		
			50-200		•			•		
	100-1000			•			•			
	Research	单通道	0.1-2.5	•			•			
			0.5-10	•			•			
			2-20		•			•		
			10-100		•			•		
			20-200		•			•		
		100-1000			•			•		
		多通道	0.5-10	•				•		
			5-100		•			•		
			10-100		•			•		
			50-100		•			•		
	20-300			•			•			
	Research Pro	单通道	20-300		•			•		
			50-1000			•			•	
		多通道	0.5-10	•				•		
			20-300	•				•		
	Research Plus	单/多通道	100-1000			•			•	
	Thermo Scientific	Finnpipette	单通道	0.2-2	•			•		
				0.3-3	•			•		
0.2-10				•			•			
0.5-5				•			•			
0.5-10				•			•			
1-10				•			•			
5-40					•			•		
5-50					•			•		
10-100					•			•		
40-200					•			•		
30-300					•			•		
100-1000						•			•	
200-1000				•			•			
多通道		1-10	•				•			
		5-50		•			•			
		50-300		•			•			
Finnpipette Digital		单通道	40-200		•					
Finnpipette electronic		多通道	1-10	•			•			
			30-300		•			•		
	100-1000				•			•		
Pipetman Classic	单通道	P2,0,1-2	•			•				
		P2,0,2-2	•			•				
		P10,0,2-10	•			•				
		P20,2-20		•			•			
		P100,10-100		•			•			
		P200,20-200		•			•			

Gilson Pipetman	Pipetman Neo	单通道	P1000,100-1000			•			•
			P2N,0.2-2	•			•		
			P10N,1-10	•			•		
			P20N,2-20		•			•	
			P100N,10-100		•			•	
			P100N,20-100		•			•	
			P200N,20-200		•			•	
	多通道	P8×20N,2-20		•			•		
	Pipetman Ultra	单通道	U200,20-200		•			•	
		多通道	P200-M8,20-200		•			•	

移液器品牌	移液器型号	通道	通道规格 (μL)	311001	312001/ 316001/310012	313001	311012	316012/ 312012	313012	314016	310012	313006/31 3016		
Axygen	Axyjet	单通道	0.1-2	•			•			•				
			0.5-10	•			•			•				
			2-20		•				•			•		
			5-50		•				•					
			10-100		•				•					
			20-200		•				•					
		100-1000			•				•				•	
		多通道	1-10	•				•			•			
			5-50		•				•					
			20-200		•				•					
Biohit	Proline	单通道	0.1-2.5	•			•			•				
			0.2-10	•			•			•				
			0.5-10	•			•			•				
			5-50		•			•						
			5-100		•			•						
			10-250		•			•						
			10-500			•				•				
			20-200		•				•					
			50-1000			•				•				
			100-1000			•				•				•
	多通道	0.2-10	•				•				•			
		5-100		•				•						
		25-250		•				•						
		50-250		•				•						
		5-100		•				•						
	e-Line	单通道	0.1-10	•				•			•			
			0.2-10	•				•			•			
			5-120		•				•					
			10-300		•				•					
			50-1000			•				•			•	
		多通道	0.2-10	•				•				•		
5-120				•				•						
10-300				•				•						
0.1-3			•				•			•				
0.5-10			•				•			•				
m-Line	单通道	10-100		•				•						
		20-100		•				•						
		100-1000			•				•			•		
	多通道	0.5-10	•				•			•				
		5-100		•				•						
		30-300		•				•						
BrandTech	Transferpet te	单通道	0.1-1	•			•			•				
			0.5-10	•				•			•			
			2-20		•				•			•		
			5-50		•				•					
			10-100		•				•					
			20-200		•				•					
			25-250		•				•					
		100-1000			•				•				•	
		多通道	0.5-10	•				•			•			
			2-20		•				•			•		
			2.5-25		•				•				•	
			5-50		•				•					
			10-100		•				•					
20-200			•				•							
Capp Aero	单/多通道	0.2-2	•				•			•				
		0.5-10	•				•			•				
		5-50		•				•						
		10-100		•				•						
		25-200		•				•						
		50-300		•				•						
		100-1000			•				•				•	
Labpette Discovery or Rinette™	单通道	0.1-2	•				•			•				
		0.5-10	•				•			•				
		2-20		•				•			•			
		10-100		•				•				•		
		20-200		•				•						
		100-1000			•				•				•	

Labnet	BioPette Plus	多通道	1-10	•			•			•		
			5-50		•			•				
			20-200		•			•				
			50-300		•			•				
	BioPette E	单通道	0.5-10	•			•			•		
			2-20	•			•			•	•	
			10-200		•			•				
		100-1000				•			•			•
		多通道	0.5-10	•			•			•		
			2-20	•			•			•	•	
10-200			•			•						
			100-1200			•			•		•	
Matrix	Impact/ Impact ²	单/多通道	2-125		•			•				
Nichryo/ Oxford	Benchmate	单通道	0.5-10	•			•			•		
			2-20		•			•			•	
			10-50		•			•				
Rainin	EDP3	单通道	E-MIC-10,0.5-10	•			•			•		
			E2-25,2.5-25		•			•				
			E2-100,10-100		•			•				
			E2-250,25-250		•			•				
			E2-1000,100-1000				•			•		
	多通道	EP-M8-10,0.5-10	•			•			•			
		EP-M8-250,10-250		•			•					
	Latch Mode	单通道	R2,0.1-2	•			•			•		
			R10,0.2-10	•			•			•		
			R20,2-20		•			•				•
R100,10-100				•			•					
R200,20-200				•			•					
R1000,100-1000				•			•			•		
Socorex	Callbra/ Accura	单通道	0.1-2	•			•			•		
			0.2-2	•			•			•		
			0.5-10	•			•			•		
			1-10	•			•			•		
			2-20		•			•				•
			10-100		•			•				
			20-200		•			•				
		100-1000				•			•			•
		多通道	1-10	•			•			•		
			5-50		•			•				
			10-100		•			•				
			20-200		•			•				
			50-200		•			•				
VWR	VWR and VWR Signature	单通道	0.1-2	•			•			•		
			0.5-10	•			•			•		
			2-20		•			•				•
			5-50		•			•				
			10-100		•			•				
			20-200		•			•				
			50-250		•			•				
		100-1000				•			•			•
		多通道	1-10	•			•			•		
			5-50		•			•				
20-200			•			•						
			50-300		•			•				
Eppendorf	Eletronic	单/多通道	0.5-10	•			•			•		
			5-100		•			•				
			2-200		•			•				
				50-1000			•			•		
	Reference	单/多通道	0.1-2.5	•			•			•		
			0.1-10	•			•			•		
			0.5-10	•			•			•		
			2-20		•			•				•
			10-100		•			•				
			20-200		•			•				
50-200		•			•							
100-1000				•			•			•		
Research	单通道	0.1-2.5	•			•			•			
		0.5-10	•			•			•			
		2-20		•			•				•	
		10-100		•			•					
	20-200		•			•						
	100-1000				•			•			•	
Research	单/多通道	0.5-10	•			•			•			
		5-100		•			•					
		10-100		•			•					

	多通道	50-100			•			•						
		20-300			•			•						
		30-300			•			•						
		50-1200				•			•			•		
	Research Pro	单通道	20-300			•			•					
		50-1000					•			•				
	多通道	0.5-10	•				•			•				
20-300		•					•			•				
Research Plus	单/多通道	100-1000				•			•			•		
Thermo Scientific	Finnpipette	单通道	0.2-2	•				•			•			
			0.3-3	•					•			•		
			0.2-10	•						•			•	
			0.5-5	•						•			•	
			0.5-10	•						•			•	
			1-10	•						•			•	
		5-40			•					•				
		5-50			•						•			
		10-100			•						•			
		40-200			•						•			
		30-300			•						•			
		100-1000					•				•		•	
	200-1000					•				•		•		
	多通道	1-10	•								•			
		5-50			•					•				
		50-300			•					•				
Finnpipette Digital	单通道	40-200			•									
Finnpipette electronic	多通道	1-10	•					•			•			
		30-300			•				•			•		
		100-1000					•			•		•		
Gilson Pipetman	Pipetman Classic	单通道	P2,0,1-2	•					•			•		
			P2,0,2-2	•						•			•	
			P10,0,2-10	•							•			•
			P20,2-20			•						•		•
			P100,10-100			•						•		
			P200,20-200			•						•		
	Pipetman Neo	单通道	P2N,0,2-2	•								•		
			P10N,1-10	•								•		
			P20N,2-20			•							•	
			P100N,10-100			•								
			P100N,20-100			•								
		P200N,20-200			•									
		多通道	P8×20N,2-20			•						•		
	Pipetman Ultra	单通道	U200,20-200			•						•		
		多通道	P200-M8,20-200			•						•		